



# Quantum Design Educational Module

## Observing the Hall Effect, Part I: Copper (Electrical Transport Option)

<http://education.qdusa.com/experiments.html>

Prof. Richard Averitt  
University of California, San Diego

*Description: The objective of this educational module is to measure the Hall voltage  $V_H$  to determine the Hall coefficient  $R_H$  of Cu, a monovalent metal.  $V_H$  in metals is typically quite small (~microvolts for reasonable values of the applied current and magnetic field), but is easily measured using VersaLab/ETO.*

### **Background:**

The Hall effect (HE)<sup>1</sup>, arises from the sideways deflection of charges in a conductor upon application of a magnetic field. This effect was discovered in 1879 by Edwin Hall and, as described below, provides a method to determine the concentration ( $n$ ) and sign (e.g. electrons or holes) of the charge carriers. The HE is now used extensively in characterizing metals, semiconductors, and ferromagnets<sup>2</sup>. In two-dimensional conductors the integer quantum Hall effect (IQHE) and fractional-QHE (FQHE) have been discovered leading to Nobel prizes in 1985 and 1998. Recently, the QHE has been observed in graphene at room temperature<sup>3</sup>. Finally, we mention that  $V_H$  is sufficiently large in semiconductors enabling the development of Hall probes, which are used as compact inexpensive magnetometers.

For a simple understanding of the HE, it is sufficient to consider the Lorentz force,


$$\mathbf{F} = e\mathbf{E} + e(\mathbf{v}_d \times \mathbf{B}) \quad (1)$$

and the current density,

$$\mathbf{J} = \frac{I}{wt} = nev_d \quad (2)$$

where (dropping the vector notation),  $E$  and  $B$  are the electric and magnetic field,  $v_d$  is the drift velocity of the charge carriers,  $e$  is the magnitude of the charge,

$J$  is the current density,  $I$  is the current,  $w$  and  $t$  are the sample width and thickness (see. Fig. 1) and  $n$  is the carrier concentration.



**Figure 1:** diagram of the Hall Effect

Figure 1 depicts a schematic of the HE for the case that the carriers are electrons (i.e. negative charges). With the magnetic field along  $z$ , the current in the negative  $x$  direction, the electron will be deflected in the positive  $y$  direction. This leads to an excess of negative charge along the top portion of the sample, with a corresponding positive density along the bottom as depicted in the Fig. 1. This charge build up leads to an electric field along  $y$  that balances the magnetic force. This results in the development of a potential difference that is transverse to the current. This transverse potential difference can be measured and is  $V_H$ , the Hall voltage.

It is easy to show using Equations 1 and 2 that,

$$V_H = wJB/ne = IB/nte = R_H IB/t \quad (3)$$

In this equation it is assumed (as in Fig. 1) that the width of the conductor is the same as the contact spacing to measure  $V_H$ . From this equation, several things are evident.  $V_H$  is linear in the current and magnetic field and inversely proportional to the carrier density and thickness. Thus, it can be expected that thinner films will yield a larger  $V_H$ , as will samples with a smaller carrier density (e.g. semiconductors).

An important quantity is the Hall coefficient  $R_H$ . For completeness, we note that

$$R_H = V_H t / IB = 1/ne \quad (4)$$

From Eqn. 4, we see that  $R_H$  can be determined exclusively from experimentally measured quantities, and that it depends on both the carrier density and sign. This equation encapsulates the power of Hall measurements. When electrons are the charge carriers,  $R_H$  is negative and when holes are the charge carriers,  $R_H$  is positive.

It is important to establish a clear convention in determining the sign of  $R_H$  since, otherwise, it is possible to misidentify the carrier type in a Hall transport measurement. Figure 1 includes labels that are consistent with the ETO transport pucks (e.g.  $I_+$ ,  $I_-$ ,  $V_+$ ,  $V_-$ ). Specifically, the field points along the z-direction that, in the VersaLab, corresponds to vertical. In addition, a positive current is from  $I_+$  to  $I_-$  and  $V_+ - V_- = V_H$  which, as shown in Fig. 1, is negative for electrons.

### Notes:

1. See, for example, Chapter 6 of *Introduction to Solid State Physics*, 7<sup>th</sup> edition, C. Kittel, Wiley and Sons, New York 1996.
2. [http://en.wikipedia.org/wiki/Hall\\_effect](http://en.wikipedia.org/wiki/Hall_effect)
3. Y. Zhang, et al, Nature 438, 201 (2005).

### Student Learning Outcomes:

- Students will learn the basic aspects of sample handling and soldering.
- Students will learn the technical aspects of performing small signal transport measurements.
- Students will learn fundamental techniques in electronic transport characterization of materials.

### Safety Information:

Before attempting to perform any parts of this student experiment, please read the entire contents of: this Educational Module, the VersaLab User's Manual (1300-001), and the Electro-Transport Option Manual (1084-700), and observe all instructions, warnings and cautions. These are provided to help you understand how to safely and properly use the equipment, perform the experiments and reach the best student learning outcomes.

Quantum Design Inc. disclaims any liability for damage to the system or injury resulting from misuse, improper operation of the system and the information contained in this Educational Module.

The following Safety warnings apply to this Educational Module. We recommend that you study them carefully and discuss the details with your instructor before starting the work:

**WARNING!**

Always use **Personal Protective Equipment (PPE)** during every step of sample preparation. Failure to do so might cause bodily harm.





### TOXIC HAZARD!

Acetone is toxic if swallowed. For more information consult the Material Safety Data Sheet available on this website:

<http://www.guidechem.com/msds/110-20-3.html>



### HOT SURFACE!

Never touch the tip of the soldering iron as it is typically at 400 °C and can cause serious burns. Please read the entire contents of the User's Manual specific to the soldering iron and solder and flux cleaners you are using, and observe all instructions, warnings and cautions. Failure to properly handle hot surfaces might cause bodily harm.

General guidelines for soldering safety can be found at:

<http://safety.eng.cam.ac.uk/procedures/Soldering/soldering-safety>


### Instructions:

Before moving on to the experimental procedure for measuring the HE in copper, we point out aspects of the VersaLab and ETO with which you should be familiar. First, read over the basic operational aspects of VersaLab as found in the user's manual (part number 1300-001, B0). This includes the helium compressor and cryocooler, magnet, puck interface, diaphragm and sorption pump, and the MultiVu software for measurement control and data acquisition.

Secondly, read over the basic features of the ETO as described in the ETO user's manual (part number 1084-700, B0). Of particular importance is the ETO hardware which includes the ETO CM-H module and EM-QN remote head which are shown in Figure 2 and the transport puck which is shown in Figure 3.



**Figure 2:** ETO electronics: module and head



**Figure 3:** Electrical transport measurement puck

We will now provide details on sample preparation and, subsequently, performing the HE measurement using the VersaLab/ETO.

Several items are needed for this experiment, which includes:

- Copper foil less than 100 microns thick
- Calipers to measure Cu foil thickness
- Cotton swabs and acetone to clear Cu foil
- Soldering iron, solder, thin gauge wire
- Cigarette paper + Apiezon H grease, or Kapton tape
- Tweezers, toothpick (for spreading H grease)
- Latex or nitrile gloves for sample handling
- ETO transport puck
- Puck wiring test station and ohm meter (to test continuity of solder joints)

Prior to performing the HE measurement with the VersaLab/ETO, it is important to prepare the sample as detailed in the following steps:

- a.) From the Cu foil, cut a square piece and make sure it can lie flat within the square region of the sample puck (see Fig. 3). Note that the measurement


puck can accommodate two measurement channels, so you may also choose to use a smaller Cu piece and simultaneously measure another sample, such as the germanium sample which is described in "Observing the Hall Effect, Part II."

b.) Using the calipers, measure and record the thickness of the Cu foil (Fig. 4).




**Figure 4**

c.) As shown in Fig. 5, clean the Cu foil using tweezers to hold the sample. A cotton swab dipped in acetone will remove any residual oils or contaminants, facilitating soldering and minimizing contamination of the VersaLab sample chamber.




**Figure 5**

d.) As shown in Fig. 6 and 7, solder a short section of wire onto each corner of the Cu foil. The excess length of wire can be cut off upon soldering the sample to the transport puck (step h, below).




**Figure 6**




**Figure 7**

e.) The transport puck is metallic in order to maintain good thermal contact between the sample and cold head. However, to prevent shorting of the sample, electrical isolation is needed. This can be achieved, for example, with Kapton tape. In Figure 8, we show an alternative option where we will use cigarette paper and Apiezon H grease. Apiezon H is a silicone-free high thermal conductivity grease suitable for high vacuum applications.




**Figure 8**

f.) As shown in Figure 9, a small section of the cigarette paper should be placed on the transport puck and covered with a thin layer of H grease.




**Figure 9**

g.) In preparation for soldering the Cu foil to the transport puck, the contacts should be pre-tinned as shown in Figure 10. In the present case, we are using channel 2 for the measurements.



**Figure 10**

h.) As shown in Figure 11 and 12, the Cu foil should be placed on the transport puck, soldering each wire and then cutting off the excess.



**Figure 11**






Figure 12

- i.) Figure 13 shows the transport puck and Cu foil attached to the puck wiring test station. Importantly, notice how the wires are soldered with  $V_+$  and  $V_-$  along one diagonal and  $I_+$  and  $I_-$  along the other diagonal (labels inset), in a configuration that is consistent with the convention for the HE described in Figure 1. Be careful to provide ample room to keep the wires from touching one another and shorting.




Figure 13

- j.) As shown in Figure 14, the puck wiring test station provides a convenient means to check the continuity of the soldering prior to insertion into the VersaLab instrument. If the solder joints are good, a resistance of less than one ohm will be measured across any of the contacts.



Figure 14

- k.) At this point the transport puck/Cu foil is ready to be inserted into the VersaLab. This is accomplished using the puck insertion tool as described in the VersaLab user's manual Figure 1-2. Upon inserting the puck, it is important to make sure the tab on the transport puck is properly aligned. The tab should face towards the front of the VersaLab, and it is possible to feel a slight click. At this point gentle downward pressure will allow for appropriate seating of the sample puck into the sample chamber.
- l.) The sample chamber can be sealed by inserting the baffle set (don't forget the o-ring) and Kwik-Flange clamp.

The rest of this experiment will utilize the MultiVu software. Please see Chapter 3 of the VersaLab manual and Chapter 4 of the ETO manual for complete details.

- m.) Activate the ETO option. Under the utilities tab, select activate option. Under the "Available Options" select "Electrical Transport" and then click the activate button. This will activate the instrument and pull up the ETO console.
- n.) On the ETO console, click the "Datafile" tab and choose a sample name and location for your data to be saved. There are default settings. Click "Change" to alter the name and location of the datafile.
- o.) While not required for this measurement, it is useful to know how to pump down the sample chamber. Under the "Instrument" tab, select "Chamber". This will pull up a dialog box. Under "Control" select "Purge/Seal" to initiate putting the sample under vacuum. The left hand side of the dialog box will show the status. This will reduce the sample vacuum to a few torr. To obtain high vacuum, click the "HiVac" tab.
- p.) The next step is to write a sequence to perform the required measurement. Under the "file" tab, click new sequence. This will open a new sequence file (e.g. Sequence1.seq) and the sequence command bar. Double click "Scan Field" to open another dialog box to enter the desired parameters.

For example, you might choose an initial field of 10000 Oe and a final field of -10000 Oe in steps of 1000 Oe. The maximum rate that can be chosen is 300 Oe/sec. Once your selections are made, click "ok" to insert the command line into the sequence file.

- q.) Returning to the sequence command bar, click "Electrical Transport" and then double click "ETO Resistance". This will pull up the "AC Resistance Measurement" box. Since we have wired up the sample on channel two, select "Enable Measurement" for Channel 2. For the excitation, choose 100mA and 21.3 Hz. Select autorange, and for the measurement configuration choose a one second averaging time and 5 for the number of measurements. Upon clicking "ok", these selections will be inserted into the sequence.
- r.) This constitutes a complete, albeit short, measurement sequence to measure the HE on copper. The sequence can be saved by selecting the "file" tab and "save as" to select a file name and location.
- s.) To initiate the measurement sequence, click the "play" button near the top of the MultiVu software interface. This will start the data acquisition with the data being save under the filename you selected in step o.) above.
- t.) At the bottom of the MultiVu software, the status of the temperature, field, and vacuum can be monitored.
- u.) To monitor the data collection in real time click the "Datafile" tab on the ETO console and select "View". This will plot the data in real time. The default is to show the time stamp on the x-axis, and to plot the resistance and phase angle of both channels one and two. To change this right click on the plot and choose "data selection" to bring up a dialog box to change the view. Change the x-axis to plot field and deselect the boxes for channel one and the phase angle of channel 2. The result is that the Hall resistance will be plotted as a function of the field. Right click on the plot and choose "autoscale" to see the complete data set.
- v.) At the end of the measurement sequence, you should obtain a plot of the Hall resistance as a function of field.
- w.) To return to the zero field condition select the "Instrument" tab and choose "field." This will open a dialog box. Enter a setpoint of 0 Oe and click "set". This will ramp the field down to zero. After venting the sample (if under vacuum), the sample can be removed.

Consider repeating this experiment at different temperatures (as suggested by Question 4 below).

### **Questions / Analysis:**

Some of the following questions are specific to the data that was obtained, while others are of a more open-ended or comparative nature.

1. From your data, determine the Hall coefficient  $R_H$  for copper. Determine the carrier density of your sample. Make certain to include the error bars for your data, along with a brief description of your error analysis.
2. Compare your results to existing values in the literature (be certain to cite your sources). Discuss possible discrepancies in terms of samples, sources of error, assumptions in your analysis, and anything else you consider as a potential contributing factor.
3. Find the value of  $R_H$  for single crystal Cu & discuss the difference with your measurements on polycrystalline samples.
4. What would you expect to happen to  $R_H$  with decreasing temperature? Why? If you would like (or have not already done so), measure  $R_H$  at 50K to see if your experimental results agree with your expectations.