

## Biowaste Treatment Log - Autoclave

**PI:** \_\_\_\_\_  
**BLDG:** \_\_\_\_\_  
**ROOM:** \_\_\_\_\_

The waste in the log below was prepared according to SOP:  
 The autoclave in room \_\_\_\_\_ was used according to SOP:

Biowaste is autoclaved at a minimum of 121°C and 15psi for no less than 30 minutes.

DATE MM/DD/YYYY	PRINTED NAME	INITIALS	QTY (lbs. or L) OF WASTE	COMMENTS